

Trimethyl Phosphate compatibility with AFLAS®



Test method Soaked into the fluid at room temperature for 72 hours.
Test fluid Trimethyl Phosphate

Test piece AFLAS® 100S (standard formulation)

Formulation	AFLAS® 100S	100
	MT-Carbon(N990)	30
	TAIC*	5
	Perkadox® 14**	1
	Sodium Stearate	1

(phr)

Cure Conditions Press molded at 170C for 20min
Post cured at 200C for 4h

Properties (before test)	AFLAS® 100S
Tensile strength [MPa]	20
Tensile Elongation [%]	230
Hardness [shore-A]	72

Trimethyl Phosphate compatibility R.T. for 72 hours	AFLAS® 100S
Change of Tensile Strength [%]	-1.1
Change of Tensile Elongation [%]	7.3
Change in Hardness [points]	0
Volume change [%]	1.2

* Triallylisocyanurate

** 1,3-bis(t-butylperoxy)-diisopropylbenzene. Perkadox® is a registered trademark of Akzo Nobel Chemicals B.V.